



ieMR Advanced

## Add, modify and remove monitor values

### SAAnaesthesia

#### Quick reference guide

There are several ways of launching the **Add Monitor Value** window. This QRG outlines these methods so that data can be documented in the **Anaesthetic Record**. **Modifications** to documented **Values** and the **Removal of Values** is also outlined below.

#### Document Monitor Values


The following outlines the three ways to launch the **Add Monitor Value** window:

- Click the **Monitor** name in the **Anaesthetic Record Grid**.




Or

- Click the **Monitor** name in the **Documentation Tabs, Monitors** tab.




Or

- Click the **Monitor** timeline for the time the user would like the **Value** documented.




- Using one of the three methods outlined above click to launch the **Add Monitor Value** window:



- Complete the required fields of **Value** and **Value Time**, a **Comment** can be added if needed.
- Click **OK**.


The **Monitor Value** entered now populates the **Anaesthetic Record Grid**:



#### Modify Documented Monitor Values

If the documented **Monitor Value** requires modification, the following steps outline how to modify documented **Monitor Values** in the **Anaesthetic Record**:

- Right click the recorded value on the timeline that requires modification, select **Modify Value**. Note: there is the option to select **Remove Value**.



The **Modify Monitor Value** window displays with the documented **Value** and **Value Time** populated:





2. Modify the data as required such as the Value, Comments section or Value Time.
3. Click OK to action the modification(s) to the Anaesthetic Record Grid. Note, there is an option to select the Remove button.


The modification populates the Anaesthetic Record Grid:



### Remove Documented Monitor Values


If the documented Monitor Value requires deletion, the following steps outline how to Remove documented Monitor Values in the Anaesthetic Record:

1. Right click the recorded value on the timeline that needs to be removed.



2. From the drop down menu displayed click to select Remove Value.

The SurgiNet: Anaesthesia window displays to confirm the deletion action:



3. Click Yes to proceed with removing the documented value. Note: clicking No returns the user to the Record.

The Monitor Value has now been removed from the Anaesthetic Record Grid:



**Note:** If a Monitor has not been Executed to the Record, the Monitor and associated data will not appear in the Anaesthetic Record Grid.

The Monitors tab of the Documentation Tabs will not display the Monitor and therefore Values cannot be recorded until the Monitor is Executed to the Anaesthetic Record.

**Note:** If a Monitor has been added to the Record and all of the Monitor Values documented are removed, there is the potential that the Monitor itself is removed from the Anaesthetic Record Grid.

In addition, the action of removing the Value(s) documented, removes the ability for the user to Add Monitor Value from the Documentation Tabs, Monitors list.

