


Interactive view: customising and configuring

Care Delivery

Quick reference guide


In **Interactive View**, **PowerChart** displays default bands, sections and fields.



It is possible to customise **Interactive View** to include additional bands, sections and fields.

Modifying the Font Size in Interactive View

1. Navigate to **Interactive View**.
2. Click the **Customise View** icon .
- A window will open to the **Customise** tab.
3. Click the **Preferences** tab.




4. Modify the font size as desired and then click **OK**.

This will apply the changes to all the **bands** in **Interactive View**.

Modifying Interactive Bands – Adding a Section

It is possible to add **Sections** into an **Interactive Band**. For example, if you wish to add the **Infant Feeding** section, you will need to customise your view.

1. Navigate to the relevant band.
2. Click on the **Customise View** icon .
3. Find the section that you want to add (this is easier to view and navigate around if you click **Collapse All**).



4. Tick the check box so that the section is on view.


Note: This section will only be saved in the current chart that you are working in, but will be viewable by all clinicians accessing that chart.

Modifying Sections – Adding a Field


It is possible to customise sections in **PowerChart** to add additional fields to input results into.

For example, **FIO2** is not a default field in the **Vital Signs** band as shown in the below screenshot, but the screen can be easily customised to add this in.

1. Click the **Customise View...** icon  - a window will open to the **Customise** tab.




2. Within the titled **On View**, items that are on view are ticked within iView. Boxes that are not ticked, such as the **FIO2**, will not currently be displayed.
3. There is also a column titled **Default Open**. Only sections that are pre-ticked will be displayed by default.



4. Tick the **On View** box beside **FIO2** to add these to the **iView** screen.

Display Name	On View	Default Open
Oxygen flow rate	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Oxygen therapy	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
FIO2	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
SBP/DBP Cuff	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

5. Click **OK**.
The selected field/s and cells will now populate on **iView**.



You can also perform a search for an item from the **Customise View** screen. Simply type what you are looking to add to the **Search for Item** field as below and a list of options will appear.


Search for Item:

In Section:


Adding an Interactive Band

1. When in a patient's chart, click on **Interactive View**.
2. On the **Toolbar**, click on **View**.
3. Click on **Layout**.
4. Click on **Navigator Bands...**






The **Navigator Documents** window will open. The list on the left hand side of the screen is the available bands to add. The list on the right is the current list visible in iView.




The up and down arrows on the right of the box can be used to change the order of the bands.

- Click **OK**.
The **Changes to Navigator Bands** window will appear.



- Select the **band** that you wish to add and then use the arrow to move it across to the **Current Document Types** list.



- Click **OK**.
You only need to close out of the patient's chart that you are currently in, not the application.

Note: This band will now be saved in your settings and will appear in all patients' charts that you open.

