

Medication/fluid view

SAAnaesthesia


Quick reference guide

The SAAnaesthesia Medication/ Fluid View is a comprehensive 'View Only' of all Medications, Intake and Output currently documented in SAAnaesthesia. The Med/Fluid View option allows the user to view total amounts and volumes charted. This window does not allow for editing or modification.

Launching Med/Fluid View

1. Click the Med/Fluid View component button in the SAAnaesthesia Toolbar.

The Medication/ Fluid View window is displayed and the screen is divided into sections - Medications, Intake and Output. If data has been documented on the Anaesthetic Record it will appear under the relevant section heading. Note, the cumulative total is in bold type.



	Start Date/Time	Stop Date/Time	Amount	Volume
Medications				
propOFol + 20mcg/ml alfentanil (IV) (No Site Specified)			160.8974 mg+	16.73 mL+
propOFol (IV) (No Site Specified)	20-Dec-2016 13:19	<Running>	160.8974 mg-	16.73 mL-
rocuronium (IV) (No Site Specified)	20-Dec-2016 16:24		160 mg	16 mL
fentanyl (IV) (No Site Specified)	20-Dec-2016 16:28		48 mg	4.8 mL
cephaZOLin (IV) (No Site Specified)	20-Dec-2016 16:38		100 mcg	2 mL
	20-Dec-2016 17:29		100 mg	2 mL
			2 g	2 vial(s)
Intakes				
compound sodium lactate (hartmanns) IV infusion 1,000 mL (IV) (No Site Specified)			102.08 mL+	
sodium chloride 0.9% IV Infusion 1,000 mL (IV) (No Site Specified)	20-Dec-2016 16:41	<Running>	102.08 mL-	
FFP - AN (IV) (No Site Specified)			0 mL+	
	20-Dec-2016 16:53	<Running>	375 mL-	
	Bag 1	20-Dec-2016 17:22	20-Dec-2016 17:25	125 mL
	Bag 2	20-Dec-2016 17:22	20-Dec-2016 17:25	125 mL
	Bag 3	20-Dec-2016 17:22	20-Dec-2016 17:25	125 mL
Outputs				
Estimated Blood Loss - AN (IV)			160 mL	
	20-Dec-2016 17:30		150 mL	
Urine Output - AN (IV)			250 mL	
	20-Dec-2016 17:30		200 mL	
Nasogastric fluid - AN (IV)			50 mL	
	20-Dec-2016 17:30		50 mL	

Within the sections the Start Date / Time, Stop Date / Time, Amount, and Volume are recorded if the item has this information available to display. If a fluid is currently running, <Running> will display in the Stop Date / Time field.

Note: The Medication/ Fluid View can be utilised as a comprehensive handover tool or as an overall view of the patient's Intake and Output status at any point in the case.

2. To exit out of the Med/Fluid View click the Close button.

