



ieMR Advanced

Skin Alterations and Incision/Wounds

Quick reference guide


Creating a Dynamic Group

1. Open the patient's chart and select *Interactive View* from the menu
2. Navigate to the section and create a clinically appropriate dynamic group




i If a patient has multiple **Skin Alterations** or **Incision/Wounds**, each will need a separate dynamic group.

3. Document the corresponding *properties* for the dynamic group, including its *location, laterality* and *location description* and click *OK*



Documenting Details

1. Double click on the *light grey bar* at the top of the left-most column of the newly created dynamic group



2. Complete all fields as clinically appropriate

i Data fields with the conditional logic icon (a diamond with a downward arrow) contain logic that opens additional fields depending on the initial selection.

i If further information about a field is required, click on the *blue hyperlinks* to see reference texts

3. Once the assessment is completed, click *Sign* ✓
4. Click *Refresh* (0 minutes ago)

Inactivating a Dynamic Group

If the *Skin Alteration* or *Wound/Incision* no longer exists, you will need to Inactivate the dynamic group

1. Open the patient's chart and select *Interactive View* from the menu
2. Navigate to the section, update clinical details and right click the name of the dynamic group label and select *Inactivate*

