



ieMR Advanced


NIKI pump – ordering/prescribing (subcutaneous infusion)

Medication Management


Quick reference guide

To order a NIKI pump, the base solution is ordered, and the Additives are prescribed within the *Continuous Details* tab In the Order Entry Field (OEF).

1. Navigate to the Orders screen and select +Add.
2. Enter “subcu 20” into the search field and select the Subcutaneous Infusion Pump 20mL (continuous).




3. Click Done – this will take you to the orders screen.




Details for sodium chloride 0.9% up to (continuous) 20 mL

Base Solution	Bag Volume	Rate	Infuse Over
sodium chloride 0.9% up to (continuous) 20 mL	20 mL	0.83 mL/hr	24 hour(s)
Additive	Additive Dose	Normalized Rate	Delivers
Weight	Weight Type	Result dt/tm	
25 kg	Manually Entered	16/03/2018 15:18:57 AEST	

4. Ensure the patient's most recent weight has pulled into the *Continuous Details* tab (if not, enter the weight manually).
5. Click on the Additive icon.
6. Search for the required medication, e.g. “Midazolam” and select the appropriate order with “subcut pump additive”.



7. Repeat step 6 for all required infusion additives.
8. Select Done once all additives have been selected.
9. Enter the total Additive Dose (free text) e.g. 2.5mg.



Note: You must specify the dose unit, e.g. mg.

10. Review the Rate and Bag Volume

Details for midazolam subcut pump additive 2.5 mg + morphine subcut pump additive

Base Solution	Bag Volume	Rate	Infuse Over
sodium chloride 0.9% up to (continuous) 20 mL	20 mL	0.83 mL/hr	24.1 hour(s)
additive	Additive Dose	Normalized Rate	Delivers
midazolam subcut pump additive	2.5 mg		EB
morphine subcut pump additive	5 mg		EB
metoclopramide subcut pump additive	10 mg		EB
total Bag Volume	20 mL		

11. Once all of the order information has been finalised, select Sign and enter your password
12. Click Refresh on the Orders Screen
13. Navigate to the MAR and Refresh
14. Recheck the correct Order and order details

Continuous Infusions	Pending Not given within 5 days.
midazolam subcut pump additive 2.5 mg	
morphine subcut pump additive 5 mg	
metoclopramide subcut pump additive 10 mg	
sodium chloride 0.9% up to (continuous) 20 mL	
Rate: 0.83 mL/hr, Total volume (mL): 20, 20 mL, Subcutaneous, start: 20-Mar-2018 18:28:00, infuse over 24.1 hour(s)	
Administration Information	
midazolam	
morphine	
metoclopramide	
Sodium Chloride 0.9% IV	

